

产品名称: COPS8 Rabbit Monoclonal Antibody
产品货号: AMRe01451



产品概述 (Summary)

产品名称 (Production Name)	COPS8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	COPS8
别名 (Alternative Names)	COP9 signalosome subunit 8; COP9; CSN8; SGN8
基因 ID (Gene ID)	10920
蛋白 ID (SwissProt ID)	Q99627.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 23 kDa

研究背景 (Background)

Component of the COP9 signalosome complex (CSN), a complex involved in various cellular and developmental processes.


产品名称: COPS8 Rabbit Monoclonal Antibody
产品货号: AMRe01451




研究领域 (Research Area)

Immunology


图片 (Image Data)



Immunocytochemistry analysis of COPS8 (green) in HeLa using COPS8 antibody, and DAPI (blue).



Western blot analysis of COPS8 in HeLa, 3T3, C6, rat Brain, K562 lysates using COPS8 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using COPS8 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: COPS8 Rabbit Monoclonal Antibody
产品货号: AMRe01451



注意事项 (Note)

For research use only .