

产品名称: Neurofascin Rabbit Monoclonal Antibody
产品货号: AMRe01464



产品概述 (Summary)

产品名称 (Production Name)	Neurofascin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Nfasc
别名 (Alternative Names)	Neurofascin
基因 ID (Gene ID)	269116.0
蛋白 ID (SwissProt ID)	Q810U3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 138 kDa; Observed MW: 186,155,75 kDa

研究背景 (Background)

Cell adhesion, ankyrin-binding protein which may be involved in neurite extension, axonal guidance, synaptogenesis,

产品名称: Neurofascin Rabbit Monoclonal Antibody
产品货号: AMRe01464




myelination and neuron-glial cell interactions.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Neurofascin in Human placenta lysates using Neurofascin antibody.



Western blot analysis of Neurofascin in mouse Brain, rat Brain, CSO7, HeLa lysates using Neurofascin antibody.

注意事项 (Note)

For research use only .