

产品名称: SNX5 Rabbit Monoclonal Antibody
产品货号: AMRe01465



产品概述 (Summary)

产品名称 (Production Name)	SNX5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SNX5
别名 (Alternative Names)	Sorting nexin-5
基因 ID (Gene ID)	27131
蛋白 ID (SwissProt ID)	Q9Y5X3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)

Involved in several stages of intracellular trafficking.


产品名称: SNX5 Rabbit Monoclonal Antibody
产品货号: AMRe01465




研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Immunocytochemistry analysis of SNX5 (green) in Hela using SNX5 antibody, and DAPI (blue).



Western blot analysis of SNX5 in mouse spleen, mouse brain lysates using SNX5 antibody.



Western blot analysis of SNX5 in K562, rat Brain, C6, Hela lysates using SNX5 antibody

注意事项 (Note)

产品名称: SNX5 Rabbit Monoclonal Antibody
产品货号: AMRe01465



For research use only .