

产品名称: PREB Rabbit Monoclonal Antibody
产品货号: AMRe01472



产品概述 (Summary)

产品名称 (Production Name)	PREB Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PREB
别名 (Alternative Names)	Prolactin regulatory element-binding protein; Mammalian guanine nucleotide exchange factor mSec12
基因 ID (Gene ID)	10113
蛋白 ID (SwissProt ID)	Q9HCU5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 45 kDa; Observed MW: 45 kDa

研究背景 (Background)

产品名称: PREB Rabbit Monoclonal Antibody
产品货号: AMRe01472




Guanine nucleotide exchange factor that specifically activates the small GTPase SAR1B.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of PREB in K562, C6, 3T3, HeLa lysates using PREB antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PREB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .