

产品名称: KAT13A Rabbit Monoclonal Antibody
产品货号: AMRe01474



产品概述 (Summary)

产品名称 (Production Name)	KAT13A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCOA1
别名 (Alternative Names)	bHLHe74; mNRC 1; NCoA 1; RIP160; SRC 1
基因 ID (Gene ID)	8648
蛋白 ID (SwissProt ID)	Q15788.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 157 kDa; Observed MW: 157 kDa

研究背景 (Background)

Nuclear receptor coactivator that directly binds nuclear receptors and stimulates the transcriptional activities in a hormone-

产品名称: KAT13A Rabbit Monoclonal Antibody
产品货号: AMRe01474




dependent fashion. Involved in the coactivation of different nuclear receptors, such as for steroids (PGR, GR and ER), retinoids (RXRs), thyroid hormone (TRs) and prostanoids (PPARs).


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of KAT13A / SRC1 in 293 lysates using KAT13A antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using SRC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .