

产品名称: BANF1 Rabbit Monoclonal Antibody
产品货号: AMRe01476



产品概述 (Summary)

产品名称 (Production Name)	BANF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BANF1
别名 (Alternative Names)	BAF; BANF1; BCRG1; BCRP1; NGPS
基因 ID (Gene ID)	8815
蛋白 ID (SwissProt ID)	O75531.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 10 kDa; Observed MW: 10 kDa

研究背景 (Background)

Plays fundamental roles in nuclear assembly, chromatin organization, gene expression and gonad development. May

产品名称: BANF1 Rabbit Monoclonal Antibody
产品货号: AMRe01476




potently compress chromatin structure and be involved in membrane recruitment and chromatin decondensation during nuclear assembly.

研究领域 (Research Area)

Microbiology

图片 (Image Data)



Western blot analysis of BANF1 in Jurkat lysates using BANF1 antibody.

注意事项 (Note)

For research use only .