

产品名称: Pirh2 Rabbit Monoclonal Antibody
产品货号: AMRe01477



产品概述 (Summary)

产品名称 (Production Name)	Pirh2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RCHY1
别名 (Alternative Names)	RCHY1; ARNIP; CHIMP; PIRH2; RNF199; ZNF363; RING finger and CHY zinc finger domain-containing protein 1; Androgen receptor N-terminal-interacting protein; CH-rich-interacting match with PLAG1; E3 ubiquitin-protein ligase Pirh2; RING finger
基因 ID (Gene ID)	25898
蛋白 ID (SwissProt ID)	Q96PM5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 30 kDa; Observed MW: 30 kDa

产品名称: Pirh2 Rabbit Monoclonal Antibody
产品货号: AMRe01477




研究背景 (Background)

Mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including p53/TP53, HDAC1 and CDKN1B. Preferentially acts on tetrameric p53/TP53.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Pirh2 in K562 lysates using Pirh2 antibody.

注意事项 (Note)

For research use only .