

产品名称: ZHX2 Rabbit Monoclonal Antibody
产品货号: AMRe01478



产品概述 (Summary)

产品名称 (Production Name)	ZHX2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZHX2
别名 (Alternative Names)	AFP regulator 1; Zinc finger and homeodomain protein 2
基因 ID (Gene ID)	22882
蛋白 ID (SwissProt ID)	Q9Y6X8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 92 kDa; Observed MW: 110 kDa

研究背景 (Background)

Acts as a transcriptional repressor.


产品名称: ZHX2 Rabbit Monoclonal Antibody
产品货号: AMRe01478




研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of ZHX2 in rat Brain, C6, HeLa lysates using ZHX2 antibody.



Immunocytochemistry analysis of ZHX2 (green) in MCF-7 using ZHX2 antibody, and DAPI (blue)

注意事项 (Note)

For research use only.