

产品名称: FDFT1 Rabbit Monoclonal Antibody
产品货号: AMRe01487



产品概述 (Summary)

产品名称 (Production Name)	FDFT1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FDFT1
别名 (Alternative Names)	DGPT; ERG9; FDFT1; SQS; Squalene synthase; SS
基因 ID (Gene ID)	2222
蛋白 ID (SwissProt ID)	P37268.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 48 kDa; Observed MW: 48 kDa

研究背景 (Background)

Critical branch point enzyme of isoprenoid biosynthesis that is thought to regulate the flux of isoprene intermediates

产品名称: FDFT1 Rabbit Monoclonal Antibody
产品货号: AMRe01487




through the sterol pathway.


研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of FDFT1 in K562, C6, HeLa lysates using FDFT1 antibody.



Immunocytochemistry analysis of FDFT1(green) in HeLa using FDFT1 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .