

产品名称: CD90 Rabbit Monoclonal Antibody
产品货号: AMRe01496



产品概述 (Summary)

产品名称 (Production Name)	CD90 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	THY1
别名 (Alternative Names)	THY1; Thy-1 membrane glycoprotein; CDw90; Thy-1 antigen; CD antigen CD90
基因 ID (Gene ID)	7070
蛋白 ID (SwissProt ID)	P04216.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 18 kDa; Observed MW: 22 kDa

研究背景 (Background)

Thy1 is a GPI-anchored, developmentally regulated protein involved in various signaling cascades that mediate neurite

产品名称: CD90 Rabbit Monoclonal Antibody
产品货号: AMRe01496




outgrowth, T cell activation, tumor suppression, apoptosis, and fibrosis. It is highly expressed on the surface of adult neurons and is thought to play a role in modulating adhesive and migratory events, such as neurite extension.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CD90/Thy1 in rat Brain lysates using CD90 antibody.

注意事项 (Note)

For research use only .