

产品名称: PARK7/DJ1 Rabbit Monoclonal Antibody
产品货号: AMRe01508



产品概述 (Summary)

产品名称 (Production Name)	PARK7/DJ1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant Rabbit Monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC-F,IHC-P,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal Antibody
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
纯化方式 (Purification)	Affinity Purified

免疫原信息 (Immunogen)

基因名 (Gene Name)	PARK7
别名 (Alternative Names)	PARK7; Protein DJ-1; Oncogene DJ1; Parkinson disease protein 7
基因 ID (Gene ID)	11315
蛋白 ID (SwissProt ID)	Q99497

产品应用 (Application)

稀释比 (Dilution Ratio)	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20
蛋白分子量 (Molecular Weight)	Calculated MW: 20 kDa; Observed MW: 20 kDa

研究背景 (Background)

Plays a role in regulating expression or stability of the mitochondrial uncoupling proteins SLC25A14 and SLC25A27 in

产品名称: PARK7/DJ1 Rabbit Monoclonal Antibody
产品货号: AMRe01508




dopaminergic neurons of the substantia nigra pars compacta and attenuates the oxidative stress induced by calcium entry into the neurons via L-type channels during pacemaking. It cooperates with Ras to increase cell transformation, it positively regulates transcription of the androgen receptor, and it may function as an indicator of oxidative stress.


研究领域 (Research Area)

Neuroscience


图片 (Image Data)



Western blot analysis of PARK7/DJ1 in K562, rat Brain, C6, 3T3, HeLa lysates using PARK7/DJ1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PARK7/DJ1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



产品名称: PARK7/DJ1 Rabbit Monoclonal Antibody
产品货号: AMRe01508



Immunocytochemistry analysis of PARK7/DJ1 (green) in Hela using PARK7/DJ1 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .