

产品名称: SOX10 Rabbit Monoclonal Antibody
产品货号: AMRe01513



产品概述 (Summary)

产品名称 (Production Name)	SOX10 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SOX10
别名 (Alternative Names)	Transcription factor SOX-10
基因 ID (Gene ID)	6663
蛋白 ID (SwissProt ID)	P56693.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa; Observed MW: 62 kDa

研究背景 (Background)

Sox10 is a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of

产品名称: SOX10 Rabbit Monoclonal Antibody
产品货号: AMRe01513




embryonic development and the determination of cell fate. Sox10 is an important regulator of neural crest and peripheral nervous system development.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of SOX10 in rat Brain, C6, HeLa lysates using SOX10 antibody.

注意事项 (Note)

For research use only .