

产品名称: Integrin alpha IIb Rabbit Monoclonal Antibody
产品货号: AMRe01525



产品概述 (Summary)

产品名称 (Production Name)	Integrin alpha IIb Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ITGA2B
别名 (Alternative Names)	ITGA2B; GP2B; ITGAB; Integrin alpha-IIb; GPalpha IIb; GPIIb; Platelet membrane glycoprotein IIb; CD41
基因 ID (Gene ID)	3674
蛋白 ID (SwissProt ID)	P08514.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 113 kDa; Observed MW: 137 kDa

研究背景 (Background)


产品名称: Integrin alpha IIb Rabbit Monoclonal Antibody
产品货号: AMRe01525

Integrin alpha-IIb/beta-3 is a receptor for fibronectin, fibrinogen, plasminogen, prothrombin, thrombospondin and vitronectin. It recognizes the sequence R-G-D in a wide array of ligands. It recognizes the sequence H-H-L-G-G-A-K-Q-A-G-D-V in fibrinogen gamma chain.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of CD41 in C6, HeLa lysates using Integrin alpha IIb antibody.

注意事项 (Note)

For research use only .