

产品名称: Cytokeratin 18 Rabbit Monoclonal Antibody
产品货号: AMRe01527



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 18 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KRT18 Cell proliferation inducing gene 46 protein; Cell proliferation inducing protein 46; Cell proliferation-inducing gene 46 protein; CK 18; CK-18; CK18; CYK 18;
别名 (Alternative Names)	CYK18; Cytokeratin 18; Cytokeratin endo B; Cytokeratin-18; K 18; K18; K1C18_HUMAN; Keratin 18; Keratin D; keratin; type I cytoskeletal 18; Keratin-18; KRT18.
基因 ID (Gene ID)	3875
蛋白 ID (SwissProt ID)	P05783.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
----------------------	--

产品名称: Cytokeratin 18 Rabbit Monoclonal Antibody
产品货号: AMRe01527



蛋白分子量 (Molecular Weight) Calculated MW: 48 kDa; Observed MW: 48 kDa


研究背景 (Background)

K18 is a type I cytoskeletal keratin. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. Keratin 18 and its filament partner keratin 8 are perhaps the most commonly found members of the intermediate filament gene family.

研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of Keratin 18 in HeLa lysates using Cytokeratin 18 antibody.

注意事项 (Note)

For research use only .