

产品名称: Rb Rabbit Monoclonal Antibody
产品货号: AMRe01530



产品概述 (Summary)

产品名称 (Production Name)	Rb Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Rb1
别名 (Alternative Names)	p110-RB1; pRb; pp105
基因 ID (Gene ID)	19645.0
蛋白 ID (SwissProt ID)	P13405.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 105 kDa; Observed MW: 105 kDa

研究背景 (Background)

Tumor suppressor that is a key regulator of the G1/S transition of the cell cycle.


产品名称: Rb Rabbit Monoclonal Antibody
产品货号: AMRe01530



研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Rb in C6, Jurkat, HeLa lysates using Rb antibody.

注意事项 (Note)

For research use only .