

产品名称: p16 ARC Rabbit Monoclonal Antibody
产品货号: AMRe01532



产品概述 (Summary)

产品名称 (Production Name)	p16 ARC Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARPC5
别名 (Alternative Names)	p16 Arc; ARPC5; ARC16; Arp2/3 protein complex subunit p16
基因 ID (Gene ID)	10092
蛋白 ID (SwissProt ID)	O15511.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 16 kDa; Observed MW: 16 kDa

研究背景 (Background)

Functions as component of the Arp2/3 complex which is involved in regulation of actin polymerization and together with

产品名称: p16 ARC Rabbit Monoclonal Antibody
产品货号: AMRe01532




an activating nucleation-promoting factor (NPF) mediates the formation of branched actin networks.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of p16 ARC in rat Brain, C6, 3T3, HeLa lysates using p16 ARC antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using p16 ARC antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .