

产品名称: Phospho-RSK4 (Ser232) Rabbit Monoclonal Antibody
产品货号: AMRe01538



产品概述 (Summary)

产品名称 (Production Name)	Phospho-RSK4 (Ser232) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPS6KA6
别名 (Alternative Names)	RPS6KA6; RSK4; Ribosomal protein S6 kinase alpha-6; S6K-alpha-6; 90 kDa ribosomal protein S6 kinase 6; p90-RSK 6; p90RSK6; Ribosomal S6 kinase 4; RSK-4; pp90RSK4
基因 ID (Gene ID)	27330
蛋白 ID (SwissProt ID)	Q9UK32.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 84 kDa; Observed MW: 84 kDa

产品名称: Phospho-RSK4 (Ser232) Rabbit Monoclonal Antibody
产品货号: AMRe01538




研究背景 (Background)

Serine/threonine kinase that may play a role in mediating the growth-factor and stress induced activation of the transcription factor CREB.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Immunocytochemistry analysis of Phospho-RSK4 (Ser232) (green) in LNcap using Phospho-RSK4 (Ser232) antibody, and DAPI(blue).



Western blot analysis of Phospho-RSK4 (Ser232) in HeLa lysates using Phospho-RSK4 (Ser232) antibody.

注意事项 (Note)

For research use only .