

产品名称: CD11b Rabbit Monoclonal Antibody
产品货号: AMRe01561



产品概述 (Summary)

产品名称 (Production Name)	CD11b Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Itgam
别名 (Alternative Names)	CD11 antigen-like family member B; CR-3 alpha chain; Cell surface glycoprotein MAC-1 subunit alpha; CD11b; Leukocyte adhesion receptor MO1
基因 ID (Gene ID)	16409.0
蛋白 ID (SwissProt ID)	P05555.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 127 kDa; Observed MW: 170 kDa

研究背景 (Background)

产品名称: CD11b Rabbit Monoclonal Antibody
产品货号: AMRe01561




Integrin ITGAM/ITGB2 is implicated in various adhesive interactions of monocytes, macrophages and granulocytes as well as in mediating the uptake of complement-coated particles and pathogens.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of CD11b in Raw264.7 lysates using CD11b antibody.

注意事项 (Note)

For research use only .