

产品名称: BNIP3 Rabbit Monoclonal Antibody
产品货号: AMRe01564



产品概述 (Summary)

产品名称 (Production Name)	BNIP3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BNIP3
别名 (Alternative Names)	BNIP3; NIP3
基因 ID (Gene ID)	664
蛋白 ID (SwissProt ID)	Q12983.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 22 kDa; Observed MW: 28 kDa

研究背景 (Background)

Apoptosis-inducing protein that, which can overcome BCL2 suppression. May play a role in repartitioning calcium between

产品名称: BNIP3 Rabbit Monoclonal Antibody
产品货号: AMRe01564




the two major intracellular calcium stores in association with BCL2.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of BNIP3 in SH-SY5Y lysates using BNIP3 antibody.



Western blot analysis of BNIP3 in rat Brain lysates using BNIP3 antibody.



Immunocytochemistry analysis of BNIP3(green) in HeLa using BNIP3 antibody, and DAPI(blue)

产品名称: BNIP3 Rabbit Monoclonal Antibody
产品货号: AMRe01564



注意事项 (Note)

For research use only .