

**产品名称: DDIT3 Rabbit Monoclonal Antibody**  
**产品货号: AMRe01568**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | DDIT3 Rabbit Monoclonal Antibody       |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC                                 |
| 种属反应性 (Reactivity)     | Human,Mouse                            |

## 产品性能 (Performance)

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                         |
| 修饰 (Modification)   | Unmodified                                                                                           |
| 同种型 (Isotype)       | IgG                                                                                                  |
| 克隆 (Clonality)      | Monoclonal                                                                                           |
| 形式 (Form)           | Liquid                                                                                               |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| 储存溶液 (Buffer)       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| 纯化方式 (Purification) | Affinity Purification                                                                                |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                                                                                        |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | DDIT3<br>DDIT3; CHOP; CHOP10; GADD153; DNA damage-inducible transcript 3 protein; DDIT-3; C/EBP-homologous protein; CHOP; C/EBP-homologous protein 10; CHOP-10; Growth arrest and DNA damage-inducible protein GADD153 |
| 别名 (Alternative Names) |                                                                                                                                                                                                                        |
| 基因 ID (Gene ID)        | 1649                                                                                                                                                                                                                   |
| 蛋白 ID (SwissProt ID)   | P35638.                                                                                                                                                                                                                |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:50-1:100             |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 19 kDa; Observed MW: 27 kDa |

产品名称: DDIT3 Rabbit Monoclonal Antibody  
产品货号: AMRe01568




## 研究背景 (Background)

Inhibits the DNA-binding activity of C/EBP and LAP by forming heterodimers that cannot bind DNA.


## 研究领域 (Research Area)

Epigenetics and Nuclear Signaling


## 图片 (Image Data)



Western blot analysis of DDIT3 in 3T3 lysates using DDIT3 antibody.



Western blot analysis of DDIT3 in K562, CHO-K1, HeLa lysates using DDIT3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using DDIT3 antibody. High-pressure and

**产品名称: DDIT3 Rabbit Monoclonal Antibody**  
**产品货号: AMRe01568**



---

temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

**注意事项 (Note)**

For research use only .