

产品名称: CD63 Rabbit Monoclonal Antibody
产品货号: AMRe01577



产品概述 (Summary)

产品名称 (Production Name)	CD63 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cd63
别名 (Alternative Names)	CD63 antigen; CD63
基因 ID (Gene ID)	12512.0
蛋白 ID (SwissProt ID)	P41731.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 40-70 kDa

研究背景 (Background)

Functions as cell surface receptor for TIMP1 and plays a role in the activation of cellular signaling cascades. Plays a role in

产品名称: CD63 Rabbit Monoclonal Antibody
产品货号: AMRe01577




the activation of ITGB1 and integrin signaling, leading to the activation of AKT, FAK/PTK2 and MAP kinases. Promotes cell survival, reorganization of the actin cytoskeleton, cell adhesion, spreading and migration, via its role in the activation of AKT and FAK/PTK2.


研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded mouse kidney using CD63 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD63 in Raw264.7 lysates using CD63 antibody.

注意事项 (Note)

For research use only .