

产品名称: CDK5RAP3 Rabbit Monoclonal Antibody
产品货号: AMRe01580



产品概述 (Summary)

产品名称 (Production Name)	CDK5RAP3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDK5RAP3
别名 (Alternative Names)	C53; IC53; LZAP; HSF-27; MST016; PP1553; OK/SW-cl.114
基因 ID (Gene ID)	80279
蛋白 ID (SwissProt ID)	Q96JB5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 57 kDa; Observed MW: 70 kDa

研究背景 (Background)

Probable tumor suppressor initially identified as a CDK5R1 interactor controlling cell proliferation (PubMed:12054757,

产品名称: CDK5RAP3 Rabbit Monoclonal Antibody
产品货号: AMRe01580




PubMed:12737517). Negatively regulates NF-kappa-B-mediated gene transcription through the control of RELA phosphorylation (PubMed:17785205, PubMed:20228063). Also regulates mitotic G2/M transition checkpoint and mitotic G2 DNA damage checkpoint (PubMed:15790566, PubMed:19223857). Through its interaction with CDKN2A/ARF and MDM2 may induce MDM2-dependent p53/TP53 ubiquitination, stabilization and activation in the nucleus, thereby promoting G1 cell cycle arrest and inhibition of cell proliferation (PubMed:16173922). May play a role in the unfolded protein response, mediating the ufmylation of multiple proteins in response to endoplasmic reticulum stress (PubMed:23152784). May also play a role in the rupture of the nuclear envelope during apoptosis (PubMed:23478299). May regulate MAPK14 activity by regulating its dephosphorylation by PPM1D/WIP1 (PubMed:21283629).


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of CDK5RAP3 in K562, 3T3 lysates using CDK5RAP3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using CDK5RAP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .

产品名称: CDK5RAP3 Rabbit Monoclonal Antibody
产品货号: AMRe01580

