

产品名称: ATF2 Rabbit Monoclonal Antibody
产品货号: AMRe01585



产品概述 (Summary)

产品名称 (Production Name)	ATF2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Atf2
别名 (Alternative Names)	Phospho-ATF2-T69/T71
基因 ID (Gene ID)	11909.0
蛋白 ID (SwissProt ID)	P16951.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 70 kDa

研究背景 (Background)

Transcriptional activator which regulates the transcription of various genes, including those involved in anti-apoptosis, cell

产品名称: ATF2 Rabbit Monoclonal Antibody
产品货号: AMRe01585




growth, and DNA damage response. Dependent on its binding partner, binds to CRE (cAMP response element) consensus sequences (5'-TGACGTCA-3') or to AP-1 (activator protein 1) consensus sequences (5'-TGACTCA-3').


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of ATF2 in 3T3, PC-12 lysates using ATF2 antibody.



Western blot analysis of ATF2 in Jurkat lysates using ATF2 antibody.

注意事项 (Note)

For research use only .