

**产品名称: 14-3-3 theta Rabbit Monoclonal Antibody**  
**产品货号: AMRe01593**



## 产品概述 (Summary)

|                        |                                         |
|------------------------|-----------------------------------------|
| 产品名称 (Production Name) | 14-3-3 theta Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody  |
| 宿主 (Host)              | Rabbit                                  |
| 应用 (Application)       | WB                                      |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                         |

## 产品性能 (Performance)

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                         |
| 修饰 (Modification)   | Unmodified                                                                                           |
| 同种型 (Isotype)       | IgG                                                                                                  |
| 克隆 (Clonality)      | Monoclonal                                                                                           |
| 形式 (Form)           | Liquid                                                                                               |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| 储存溶液 (Buffer)       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| 纯化方式 (Purification) | Affinity Purification                                                                                |

## 免疫原信息 (Immunogen)

|                        |                                                                                     |
|------------------------|-------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | YWHAQ                                                                               |
| 别名 (Alternative Names) | YWHAQ; 14-3-3 protein theta; 14-3-3 protein T-cell; 14-3-3 protein tau; Protein HS1 |
| 基因 ID (Gene ID)        | 10971                                                                               |
| 蛋白 ID (SwissProt ID)   | P27348.                                                                             |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000                            |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 28 kDa; Observed MW: 28 kDa |

## 研究背景 (Background)

**产品名称: 14-3-3 theta Rabbit Monoclonal Antibody**  
**产品货号: AMRe01593**




Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner.


## 研究领域 (Research Area)

Neuroscience

## 图片 (Image Data)



Western blot analysis of 14-3-3 Theta in rat brain, rat lung lysates using 14-3-3 theta antibody.



Western blot analysis of 14-3-3 Theta/Tau in 3T3, HeLa lysates using 14-3-3 Theta/Tau antibody.

## 注意事项 (Note)

For research use only .