

产品名称: 15 PGDH Rabbit Monoclonal Antibody
产品货号: AMRe01594



产品概述 (Summary)

产品名称 (Production Name)	15 PGDH Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HPGD
别名 (Alternative Names)	15-PGDH; Hpgd; PGDH; PGDH1; PHOAR1; SDR36C1
基因 ID (Gene ID)	3248
蛋白 ID (SwissProt ID)	P15428.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa; Observed MW: 25 kDa

研究背景 (Background)

Prostaglandin inactivation. Contributes to the regulation of events that are under the control of prostaglandin levels.

产品名称: 15 PGDH Rabbit Monoclonal Antibody
产品货号: AMRe01594




Catalyzes the NAD-dependent dehydrogenation of lipoxin A4 to form 15-oxo-lipoxin A4. Inhibits in vivo proliferation of colon cancer cells.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of 15PGDH in SW480 lysates using 15 PGDH antibody.

注意事项 (Note)

For research use only .