

产品名称: AAMP Rabbit Monoclonal Antibody
产品货号: AMRe01597



产品概述 (Summary)

产品名称 (Production Name)	AAMP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AAMP
别名 (Alternative Names)	Angio-associated migratory cell protein
基因 ID (Gene ID)	14
蛋白 ID (SwissProt ID)	Q13685.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)

Plays a role in angiogenesis and cell migration. In smooth muscle cell migration, may act through the RhoA pathway.


产品名称: AAMP Rabbit Monoclonal Antibody
产品货号: AMRe01597




研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of AAMP in mouse heart lysates using AAMP antibody.



Western blot analysis of AAMP in C6, Hela lysates using AAMP antibody.

注意事项 (Note)

For research use only .