

产品名称: ACTN3 Rabbit Monoclonal Antibody
产品货号: AMRe01610



产品概述 (Summary)

产品名称 (Production Name)	ACTN3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACTN3
别名 (Alternative Names)	ACTN3; Alpha-actinin-3; Alpha-actinin skeletal muscle isoform 3; F-actin cross-linking protein
基因 ID (Gene ID)	89
蛋白 ID (SwissProt ID)	Q08043.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 103 kDa; Observed MW: 103 kDa

研究背景 (Background)

产品名称: ACTN3 Rabbit Monoclonal Antibody
产品货号: AMRe01610




F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of ACTN3 in rat skeletal muscle lysates using ACTN3 antibody.



Western blot analysis of ACTN3 in 3T3 lysates using ACTN3 antibody.

注意事项 (Note)

For research use only .