

产品名称: ADH5 Rabbit Monoclonal Antibody
产品货号: AMRe01611



产品概述 (Summary)

产品名称 (Production Name)	ADH5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ADH5
别名 (Alternative Names)	FDH; ADHX; ADH-3; FALDH; GSNOR; GSH-FDH; HEL-S-60p
基因 ID (Gene ID)	128
蛋白 ID (SwissProt ID)	P11766.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 40 kDa; Observed MW: 40 kDa

研究背景 (Background)

Class-III ADH is remarkably ineffective in oxidizing ethanol, but it readily catalyzes the oxidation of long-chain primary

产品名称: ADH5 Rabbit Monoclonal Antibody
产品货号: AMRe01611




alcohols and the oxidation of S-(hydroxymethyl) glutathione. Miscellaneous There are 7 different ADH's isozymes in human: three belongs to class-I: alpha, beta, and gamma, one to class-II: pi, one to class-III: chi, one to class-IV: ADH7 and one to class-V: ADH6.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of ADH5 in MCF-7, HepG2, K562 lysates using ADH5 antibody.

注意事项 (Note)

For research use only .