

产品名称: AGXT Rabbit Monoclonal Antibody
产品货号: AMRe01620



产品概述 (Summary)

产品名称 (Production Name)	AGXT Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AGXT
别名 (Alternative Names)	AGT; PH1; SPT; AGT1; SPAT; TLH6; AGXT1
基因 ID (Gene ID)	189
蛋白 ID (SwissProt ID)	P21549.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 43 kDa; Observed MW: 43 kDa

研究背景 (Background)

Peroxisomal aminotransferase that catalyzes the transamination of glyoxylate to glycine and contributes to the glyoxylate

产品名称: AGXT Rabbit Monoclonal Antibody
产品货号: AMRe01620




detoxification.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of AGXT in mouse liver, rat liver lysates using AGXT antibody.

注意事项 (Note)

For research use only .