

产品名称: ALAS2 Rabbit Monoclonal Antibody
产品货号: AMRe01623



产品概述 (Summary)

产品名称 (Production Name)	ALAS2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALAS2
别名 (Alternative Names)	ASB; ANH1; XLSA; ALASE; XLDPP; XLEPP; ALAS-E; SIDBA1
基因 ID (Gene ID)	212
蛋白 ID (SwissProt ID)	P22557.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 65 kDa; Observed MW: 65 kDa

研究背景 (Background)

MiscellaneousThere are two delta-ALA synthases in vertebrates: an erythroid- specific form and one (housekeeping) which

产品名称: ALAS2 Rabbit Monoclonal Antibody
产品货号: AMRe01623




is expressed in all tissues.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of ALAS2 in K562, C6, 3T3 lysates using ALAS2 antibody.

注意事项 (Note)

For research use only .