

产品名称: ALDH2 Rabbit Monoclonal Antibody
产品货号: AMRe01625



产品概述 (Summary)

产品名称 (Production Name)	ALDH2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALDH2
别名 (Alternative Names)	ALDH2; ALDM; Aldehyde dehydrogenase; mitochondrial; ALDH class 2; ALDH-E2; ALDHI
基因 ID (Gene ID)	217
蛋白 ID (SwissProt ID)	P05091.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 56 kDa; Observed MW: 56 kDa

研究背景 (Background)

产品名称: ALDH2 Rabbit Monoclonal Antibody
产品货号: AMRe01625




Alcohol metabolism; ethanol degradation; acetate from ethanol: step 2/2.


研究领域 (Research Area)

Tags & Cell Markers


图片 (Image Data)



Western blot analysis of ALDH2 in mouse liver lysates using ALDH2 antibody.



Western blot analysis of ALDH2 in K562 lysates using ALDH2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using ALDH2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: ALDH2 Rabbit Monoclonal Antibody
产品货号: AMRe01625



注意事项 (Note)

For research use only .