

产品名称: ALIX Rabbit Monoclonal Antibody
产品货号: AMRe01628



产品概述 (Summary)

产品名称 (Production Name)	ALIX Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Hamster,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PDCD6IP
别名 (Alternative Names)	AIP1; ALIX; HP95; DRIP4
基因 ID (Gene ID)	10015
蛋白 ID (SwissProt ID)	Q8WUM4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 96 kDa; Observed MW: 96 kDa

研究背景 (Background)

Class E VPS protein involved in concentration and sorting of cargo proteins of the multivesicular body (MVB) for

产品名称: ALIX Rabbit Monoclonal Antibody
产品货号: AMRe01628




incorporation into intraluminal vesicles (ILVs) that are generated by invagination and scission from the limiting membrane of the endosome.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of ALIX in CHO-K1, rat Brain, Jurkat lysates using ALIX antibody.

注意事项 (Note)

For research use only .