

产品名称: ALS2CR8 Rabbit Monoclonal Antibody
产品货号: AMRe01640



产品概述 (Summary)

产品名称 (Production Name)	ALS2CR8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CARF
别名 (Alternative Names)	ALS2CR8; NYD-SP24
基因 ID (Gene ID)	79800
蛋白 ID (SwissProt ID)	Q8N187.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 81 kDa; Observed MW: 81 kDa

研究背景 (Background)

Acts as a transcriptional activator that mediates the calcium- and neuron-selective induction of BDNF exon III transcription.

产品名称: ALS2CR8 Rabbit Monoclonal Antibody
产品货号: AMRe01640




Binds to the consensus calcium-response element CaRE1 5'-CTATTTCGAG-3' sequence.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of ALS2CR8 in 293T lysates using ALS2CR8 antibody.

注意事项 (Note)

For research use only .