

产品名称: Angiotensinogen Rabbit Monoclonal Antibody
产品货号: AMRe01648



产品概述 (Summary)

产品名称 (Production Name)	Angiotensinogen Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Agt
别名 (Alternative Names)	Angiotensinogen; Serpin A8; AGT; SERPINA8; Des-Asp[1]-angiotensin II; Serine (or cysteine) proteinase inhibitor
基因 ID (Gene ID)	11606.0
蛋白 ID (SwissProt ID)	P11859.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 52 kDa

研究背景 (Background)

产品名称: Angiotensinogen Rabbit Monoclonal Antibody
产品货号: AMRe01648




Essential component of the renin-angiotensin system (RAS), a potent regulator of blood pressure, body fluid and electrolyte homeostasis.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Angiotensinogen in mouse liver lysates using Angiotensinogen antibody.

注意事项 (Note)

For research use only .