

产品名称: AP2 alpha Rabbit Monoclonal Antibody
产品货号: AMRe01655



产品概述 (Summary)

产品名称 (Production Name)	AP2 alpha Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TFAP2A
别名 (Alternative Names)	AP-2; BOFS; AP2TF; TFAP2; AP-2alpha
基因 ID (Gene ID)	7020
蛋白 ID (SwissProt ID)	P05549.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 48 kDa; Observed MW: 48 kDa

研究背景 (Background)

Sequence-specific DNA-binding protein that interacts with inducible viral and cellular enhancer elements to regulate

产品名称: AP2 alpha Rabbit Monoclonal Antibody
产品货号: AMRe01655




transcription of selected genes. AP-2 factors bind to the consensus sequence 5'-GCCNNNGGC-3' and activate genes involved in a large spectrum of important biological functions including proper eye, face, body wall, limb and neural tube development.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of AP2alpha in C6 lysates using AP2 alpha antibody.

注意事项 (Note)

For research use only .