

产品名称: ApoER2 Rabbit Monoclonal Antibody
产品货号: AMRe01664



产品概述 (Summary)

产品名称 (Production Name)	ApoER2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LRP8
别名 (Alternative Names)	LRP8; APOER2; Low-density lipoprotein receptor-related protein 8; LRP-8; Apolipoprotein E receptor 2
基因 ID (Gene ID)	7804
蛋白 ID (SwissProt ID)	Q14114.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 106 kDa; Observed MW: 130,106 kDa

研究背景 (Background)

产品名称: ApoER2 Rabbit Monoclonal Antibody
产品货号: AMRe01664




Cell surface receptor for Reelin (RELN) and apolipoprotein E (apoE)-containing ligands. LRP8 participates in transmitting the extracellular Reelin signal to intracellular signaling processes, by binding to DAB1 on its cytoplasmic tail.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of ApoER2 in Human placenta lysates using ApoER2 antibody.



Western blot analysis of ApoER2 in CHO-K1, C6 lysates using ApoER2 antibody.

注意事项 (Note)

For research use only .