

产品名称: Aquaporin 5 Rabbit Monoclonal Antibody
产品货号: AMRe01670



产品概述 (Summary)

产品名称 (Production Name)	Aquaporin 5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AQP5
别名 (Alternative Names)	AQP5; Aquaporin-5; AQP-5
基因 ID (Gene ID)	362
蛋白 ID (SwissProt ID)	P55064.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa; Observed MW: 28 kDa

研究背景 (Background)

Forms a water-specific channel. Implicated in the generation of saliva, tears, and pulmonary secretions. Required for TRPV4

产品名称: Aquaporin 5 Rabbit Monoclonal Antibody
产品货号: AMRe01670




activation by hypotonicity (PubMed:16571723). Together with TRPV4, controls regulatory volume decrease in salivary epithelial cells (PubMed:16571723).

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Aquaporin 5 in MCF-7 lysates using Aquaporin 5 antibody.

注意事项 (Note)

For research use only .