

产品名称: Actin Related Protein 3 Rabbit Monoclonal Antibody
产品货号: AMRe01676



产品概述 (Summary)

产品名称 (Production Name)	Actin Related Protein 3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACTR3
别名 (Alternative Names)	ACTR3; ARP3; Actin-related protein 3; Actin-like protein 3
基因 ID (Gene ID)	10096
蛋白 ID (SwissProt ID)	P61158.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)

产品名称: Actin Related Protein 3 Rabbit Monoclonal Antibody
产品货号: AMRe01676




The Arp2/3 complex involved in Golgi polarization in NIH3T3 cells. In a different model it has been shown that the Arp2/3 complex is necessary for neutrophil chemotaxis and phagocytosis.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded mouse colon using Actin Related Protein 3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Arp3 in K562, rat Brain, C6, 3T3, Hela lysates using Arp3 antibody.



产品名称: Actin Related Protein 3 Rabbit Monoclonal Antibody
产品货号: AMRe01676



Immunocytochemistry analysis of Arp3(green) in Hela cells using Arp3 antibody,and DAPI(blue)

注意事项 (Note)

For research use only .