

产品名称: AS160 Rabbit Monoclonal Antibody
产品货号: AMRe01681



产品概述 (Summary)

产品名称 (Production Name)	AS160 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TBC1D4
别名 (Alternative Names)	AS160; NIDDM5
基因 ID (Gene ID)	9882
蛋白 ID (SwissProt ID)	O60343.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 147 kDa; Observed MW: 160 kDa

研究背景 (Background)

May act as a GTPase-activating protein for RAB2A, RAB8A, RAB10 and RAB14. Isoform 2 promotes insulin-induced glucose

产品名称: AS160 Rabbit Monoclonal Antibody
产品货号: AMRe01681




transporter SLC2A4/GLUT4 translocation at the plasma membrane, thus increasing glucose uptake.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of AS160 in 293 lysates using AS160 antibody.

注意事项 (Note)

For research use only .