

产品名称: Ataxin 1 Rabbit Monoclonal Antibody
产品货号: AMRe01686



产品概述 (Summary)

产品名称 (Production Name)	Ataxin 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATXN1
别名 (Alternative Names)	ATX1; SCA1; D6S504E
基因 ID (Gene ID)	6310
蛋白 ID (SwissProt ID)	P54253.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 87 kDa; Observed MW: 105 kDa

研究背景 (Background)

Chromatin-binding factor that repress Notch signaling in the absence of Notch intracellular domain by acting as a CBF1

产品名称: Ataxin 1 Rabbit Monoclonal Antibody
产品货号: AMRe01686




corepressor. Binds to the HEY promoter and might assist, along with NCOR2, RBPJ-mediated repression. Binds RNA in vitro. May be involved in RNA metabolism (PubMed:21475249). In concert with CIC and ATXN1L, involved in brain development .

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Ataxin 1 in 3T3 lysates using Ataxin 1 antibody.

注意事项 (Note)

For research use only .