

**产品名称: Phospho-ATP Citrate Synthase (Thr447/Ser451)
Rabbit Monoclonal Antibody**
产品货号: AMRe01694



产品概述 (Summary)

产品名称 (Production Name)	Phospho-ATP Citrate Synthase (Thr447/Ser451) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACLY
别名 (Alternative Names)	ACLY; ATP-citrate synthase; ATP-citrate; pro-S-)-lyase; ACL; Citrate cleavage enzyme
基因 ID (Gene ID)	47
蛋白 ID (SwissProt ID)	P53396.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 121 kDa; Observed MW: 121 kDa

**产品名称: Phospho-ATP Citrate Synthase (Thr447/Ser451)
Rabbit Monoclonal Antibody**
产品货号: AMRe01694




研究背景 (Background)

ATP citrate lyase is the primary enzyme responsible for the synthesis of cytosolic acetyl-CoA in many tissues. The enzyme is a tetramer (relative molecular weight approximately 440,000) of apparently identical subunits. It catalyzes the formation of acetyl-CoA and oxaloacetate from citrate and CoA with a concomitant hydrolysis of ATP to ADP and phosphate. The product, acetyl-CoA, serves several important biosynthetic pathways, including lipogenesis and cholesterologenesis.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Phospho-ATP citrate lyase (Thr447+Ser451) in K562 lysates using Phospho-ATP Citrate Synthase (Thr447/Ser451) antibody.

注意事项 (Note)

For research use only .