

产品名称: AZI2 Rabbit Monoclonal Antibody
产品货号: AMRe01706



产品概述 (Summary)

产品名称 (Production Name)	AZI2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AZI2
别名 (Alternative Names)	AZ2; NAP1; TILP
基因 ID (Gene ID)	64343
蛋白 ID (SwissProt ID)	Q9H6S1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 45 kDa; Observed MW: 45 kDa

研究背景 (Background)

Adapter protein which binds TBK1 and IKKε playing a role in antiviral innate immunity. Activates serine/threonine-protein

产品名称: AZI2 Rabbit Monoclonal Antibody
产品货号: AMRe01706




kinase TBK1 and facilitates its oligomerization. Enhances the phosphorylation of NF-kappa-B p65 subunit RELA by TBK1. Promotes TBK1-induced as well as TNF-alpha or PMA-induced activation of NF-kappa-B. Participates in IFNB promoter activation via TICAM1.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of AZI2 in K562, rat Brain, C6, HeLa lysates using AZI2 antibody.

注意事项 (Note)

For research use only .