

产品名称: BAG3 Rabbit Monoclonal Antibody
产品货号: AMRe01713



产品概述 (Summary)

产品名称 (Production Name)	BAG3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BAG3
别名 (Alternative Names)	BAG3; BIS; BAG family molecular chaperone regulator 3; BAG-3; Bcl-2-associated athanogene 3; Bcl-2-binding protein Bis; Docking protein CAIR-1
基因 ID (Gene ID)	9531
蛋白 ID (SwissProt ID)	O95817.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 62 kDa; Observed MW: 80 kDa

研究背景 (Background)

产品名称: BAG3 Rabbit Monoclonal Antibody
产品货号: AMRe01713




Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release. Has anti-apoptotic activity.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Bag3 in C6, 3T3, HeLa lysates using BAG3 antibody.

注意事项 (Note)

For research use only .