

产品名称: BAP31 Rabbit Monoclonal Antibody
产品货号: AMRe01716



产品概述 (Summary)

产品名称 (Production Name)	BAP31 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BCAP31 BCAP31; BAP31; DXS1357E; B-cell receptor-associated protein 31; BCR-associated protein 31; Bap31; 6C6-AG tumor-associated antigen; Protein CDM;
别名 (Alternative Names)	p28
基因 ID (Gene ID)	10134
蛋白 ID (SwissProt ID)	P51572.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa; Observed MW: 28 kDa

产品名称: BAP31 Rabbit Monoclonal Antibody
产品货号: AMRe01716




研究背景 (Background)

May play a role in anterograde transport of membrane proteins from the endoplasmic reticulum to the Golgi. May be involved in CASP8-mediated apoptosis.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of BAP31 in K562, HeLa lysates using BAP31 antibody.



Immunocytochemistry analysis of BAP31(green) in HeLa using BAP31 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using BAP31 antibody. High-pressure and

产品名称: BAP31 Rabbit Monoclonal Antibody
产品货号: AMRe01716



temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .