

产品名称: BNIP3L Rabbit Monoclonal Antibody
产品货号: AMRe01732



产品概述 (Summary)

产品名称 (Production Name)	BNIP3L Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BNIP3L
别名 (Alternative Names)	Adenovirus E1B19k binding protein B5; BNIP3a; BNIP3H; BNIP3L; BNIP3L protein; NIP3 like protein X; NIP3L; NIX
基因 ID (Gene ID)	665
蛋白 ID (SwissProt ID)	O60238.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa; Observed MW: 38 kDa

研究背景 (Background)

产品名称: BNIP3L Rabbit Monoclonal Antibody
产品货号: AMRe01732




Induces apoptosis. Interacts with viral and cellular anti-apoptosis proteins. Can overcome the suppressors BCL-2 and BCL-XL, although high levels of BCL-XL expression will inhibit apoptosis. May function as a tumor suppressor. Inhibits apoptosis induced by BNIP3.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of BNIP3L in Jurkat lysates using BNIP3L antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using BNIP3L antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .