

产品名称: Bora Rabbit Monoclonal Antibody
产品货号: AMRe01734



产品概述 (Summary)

产品名称 (Production Name)	Bora Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BORA
别名 (Alternative Names)	C13orf34
基因 ID (Gene ID)	79866
蛋白 ID (SwissProt ID)	Q6PGQ7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 61 kDa; Observed MW: 61 kDa

研究背景 (Background)

Required for the activation of AURKA at the onset of mitosis.


产品名称: Bora Rabbit Monoclonal Antibody
产品货号: AMRe01734




研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Bora in mouse liver lysates using Bora antibody.



Western blot analysis of Bora in K562, C6 lysates using Bora antibody.

注意事项 (Note)

For research use only .