

产品名称: C4 Binding Protein Rabbit Monoclonal Antibody
产品货号: AMRe01746



产品概述 (Summary)

产品名称 (Production Name)	C4 Binding Protein Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	C4BPB
别名 (Alternative Names)	C4BP; C4BPB; PRP
基因 ID (Gene ID)	725
蛋白 ID (SwissProt ID)	P20851.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa; Observed MW: 40 kDa

研究背景 (Background)

Controls the classical pathway of complement activation. It binds as a cofactor to C3b/C4b inactivator (C3bINA), which then

产品名称: C4 Binding Protein Rabbit Monoclonal Antibody
产品货号: AMRe01746




hydrolyzes the complement fragment C4b.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of C4 binding protein in Human serum lysates using C4 Binding Protein antibody.

注意事项 (Note)

For research use only .