

产品名称: CaMKII delta Rabbit Monoclonal Antibody
产品货号: AMRe01751



产品概述 (Summary)

产品名称 (Production Name)	CaMKII delta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CAMK2D
别名 (Alternative Names)	CAM kinase 2 delta; CAM kinase II delta; CaM kinase II delta subunit; CaM kinase II subunit delta; CAMK 2d; CaMK-II delta subunit; CAMK2D; CAMKI
基因 ID (Gene ID)	817
蛋白 ID (SwissProt ID)	Q13557.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 56 kDa; Observed MW: 56 kDa

研究背景 (Background)

产品名称: CaMKII delta Rabbit Monoclonal Antibody
产品货号: AMRe01751




CaM-kinase II (CAMK2) is a prominent kinase in the central nervous system that may function in long-term potentiation and neurotransmitter release.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of CaMKII delta in rat Brain, C6, HeLa lysates using CaMKII delta antibody.

注意事项 (Note)

For research use only .