

产品名称: Caveolin 2 Rabbit Monoclonal Antibody
产品货号: AMRe01769



产品概述 (Summary)

产品名称 (Production Name)	Caveolin 2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CAV2
别名 (Alternative Names)	CAV2; Caveolae protein 20 Kd; Caveolin 2; Caveolin 2 isoform a and b; Caveolin-2
基因 ID (Gene ID)	858
蛋白 ID (SwissProt ID)	P51636.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 18 kDa; Observed MW: 18 kDa

研究背景 (Background)

产品名称: Caveolin 2 Rabbit Monoclonal Antibody
产品货号: AMRe01769




May act as a scaffolding protein within caveolar membranes. Interacts directly with G-protein alpha subunits and can functionally regulate their activity. Acts as an accessory protein in conjunction with CAV1 in targeting to lipid rafts and driving caveolae formation.


研究领域 (Research Area)

Tags & Cell Markers


图片 (Image Data)



Immunocytochemistry analysis of Caveolin 2 (green) in A549 using Caveolin 2 antibody, and DAPI (blue).



Western blot analysis of Caveolin2 in HeLa, CHO-K1 lysates using Caveolin2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using Caveolin2 antibody. High-pressure

产品名称: Caveolin 2 Rabbit Monoclonal Antibody
产品货号: AMRe01769



and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .