

产品名称: CD168 Rabbit Monoclonal Antibody
产品货号: AMRe01778



产品概述 (Summary)

产品名称 (Production Name)	CD168 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HMMR HMMR; IHABP; RHAMM; Hyaluronan mediated motility receptor; Intracellular
别名 (Alternative Names)	hyaluronic acid-binding protein; Receptor for hyaluronan-mediated motility; CD antigen CD168
基因 ID (Gene ID)	3161
蛋白 ID (SwissProt ID)	O75330.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 84 kDa; Observed MW: 84 kDa

产品名称: CD168 Rabbit Monoclonal Antibody
产品货号: AMRe01778




研究背景 (Background)

Involved in cell motility. When hyaluronan binds to HMMR, the phosphorylation of a number of proteins, including PTK2/FAK1 occurs. May also be involved in cellular transformation and metastasis formation, and in regulating extracellular-regulated kinase (ERK) activity.

研究领域 (Research Area)

Stem Cells

图片 (Image Data)



Western blot analysis of CD168 in K562 lysates using CD168 antibody.

注意事项 (Note)

For research use only .